

Artemis DLL Version History

Version: 708

Created: 27th September 2016

By: Chris Golden

Change	Description
Added Camera	Added 'Synoptics USB 2.0 12.0' camera to list

Version: 709

Created: 13th October 2016

By: Chris Golden

Change	Description
IC24 Camera Flush	If connecting to IC24 fails, we attempt to flush all the data from the camera before retrying
Infinity EEPROM Read / Write	Added ability to read and write to eeprom of infinity

Version: 710

Created: 13th October 2016

By: Chris Golden

Change	Description
GP Camera	Abort Exposure now works properly

Version: 711

Created: 21st October 2016

By: Chris Golden

Change	Description
Infinity EEprom Read / Write	Can only happen on local connections

Version: 712

Created: 28th October 2016

By: Chris Golden

Change	Description
460ex Colour offsets	Now sets bayer offsets properly

Version: 713

Created: 1st November 2016

By: Chris Golden

Change	Description
AtikMessagePing	Now set to 0 bytes.
Artemis Temperature Messages	Now only sent if we have received data within the last 10 seconds. This is done to prevent the send buffer from filling up when the client is idle.

Version: 714

Created: 16th November 2016

By: Chris Golden

Change	Description
SonySci Bin Bug	Max Y bin now set correctly

Version: 715

Created: 7th December 2016

By: Chris Golden

Change	Description
HSC Window Heater	Added HSC window heater code
CameraStub Bug	CameraStubs no longer cleared in device refresh

Version: 716

Created: 14th December 2016

By: Chris Golden

Change	Description
Temperature Bug Fix	Temperature now send to client properly
Window Heater Bug	Values greater than 128 are no longer negative.

Version: 717

Created: 21st December 2016

By: Chris Golden

Change	Description
SciCam Improvement	SetGain method now uses sanity checking to prevent errors.

Version: 718

Created: 17th January 2017

By: Chris Golden

Change	Description
VID and PID	Can now set VID and PID of cameras

Version: 719

Created: 23rd January 2017

By: Chris Golden

Change	Description
Sony Sci Cooling Power	Can now use power based cooling on sony sci cams.

Version: 720

Created: 16th Feb 2017

By: Chris Golden

Change	Description
16200 Column Correct	Can now set bad columns on 16200 cameras
CMOS Panasonic	Added to dll (not quite finished)

Version: 721

Created: 21st Feb 2017

By: Chris Golden

Change	Description
Exposure Thread Clash	Same serial numbers no longer causes clash
Disconnect Camera Crash	Service no longer crashes when camera disconnect multiple times

Version: 722

Created: 23rd Feb 2017

By: Chris Golden

Change	Description
Firmware Uploader	IC24 Cameras
ArtemisColumnRepairCanFixColumns	Added to make things easier